



Introducing our new range of centre folding rollators.

Features include:

- ✓ Compact X folding mechanism for easy storage.
- ✓ Stands independently when folded.
- ✓ Perfectly balanced frame when lifting.
- ✓ Effortless braking system.
- ✓ Simple handle height adjustment.
- ✓ Basket & shopping bag available.
- ✓ Loading capacity : 113kg / 250lbs
- ✓ Available in 3 seat heights:
 - * STANDARD (61cm / 24")
 - * LOW (54cm / 21.3")
 - * SUPER LOW (48cm / 18.9")

Exclusively for you

ROLLATOR

X Folding - Simple handle - Light weight

RA68



RG77



RA66



Trolli Master Sdn. Bhd (209382-P)
 Lot 525, Mukim Telok Panglima Garang,
 42500 Daerah Kuala Langat,
 Selangor Darul Ehsan

Tel : +60(3) 3191 4077, 3191 2837
 Fax : +60(3) 3191 2152
 Email : tm@trollimaster.com.my
 Website : www.trollimaster.com



RA60

RA60	Handle height (min - max)	Seat height	Length	Width	Width between handle	Net weight	Weight capacity
STANDARD	81 - 101 cm / 32" - 39.8"	61cm / 24"	65cm / 25.6"	65cm / 25.6"	48cm / 18.9"	7.2kg / 15.8lbs	113kg / 250lbs
LOW	74 - 94cm / 29" - 37"	54cm / 21.3"	65cm / 25.6"	65cm / 25.6"	48cm / 18.9"	7.1kg / 15.6lbs	113kg / 250lbs
SUPER LOW	70 - 90cm / 27.6" - 35.4"	48cm / 18.9"	65cm / 25.6"	65cm / 25.6"	48cm / 18.9"	7.0kg / 15.4lbs	113kg / 250lbs

- 7 handle height adjustments.
- Available with basket or shopping bag.
- Maximum basket and shopping bag load 5kg / 11lbs.
- 200mm EVA wheel.

Highlight Features



Large padded seat



Simple handle adjustment



Cable hidden braking system

NEW



RA66

RA66	Handle height (min - max)	Seat height	Length	Width	Width between handle	Net weight	Weight capacity
STANDARD	81 - 101 cm / 32" - 39.8"	61cm / 24"	65cm / 25.6"	65cm / 25.6"	48cm / 18.9"	6.2kg / 13.4lbs	113kg / 250lbs
LOW	74 - 94cm / 29" - 37"	54cm / 21.3"	65cm / 25.6"	65cm / 25.6"	48cm / 18.9"	6.1kg / 13.2lbs	113kg / 250lbs
SUPER LOW	70 - 90cm / 27.6" - 35.4"	48cm / 18.9"	65cm / 25.6"	65cm / 25.6"	48cm / 18.9"	6.0kg / 13.0lbs	113kg / 250lbs

- Cable hidden braking system.
- 7 handle height adjustments.
- Available with shopping bag only.
- Maximum shopping bag load 5kg / 11lbs.
- 200mm EVA wheel.

Highlight Features



Net seat



Secondary frame lock



Folds compactly for storage. 25cm depth when folded

NEW



RA68

RA68	Handle height (min - max)	Seat height	Length	Width	Width between handle	Net weight	Weight capacity
STANDARD	79 - 99cm / 31" - 39"	58cm / 23"	65cm / 25.6"	65cm / 25.6"	48cm / 18.9"	6.9kg / 15.8lbs	113kg / 250lbs
LOW	71 - 91cm / 28" - 36"	51cm / 20"	65cm / 25.6"	65cm / 25.6"	48cm / 18.9"	6.8kg / 15.6lbs	113kg / 250lbs
SUPER LOW	68.5 - 87.5cm / 27" - 34.5"	46cm / 18"	65cm / 25.6"	65cm / 25.6"	48cm / 18.9"	6.7kg / 15.4lbs	113kg / 250lbs

- 7 handle height adjustment
- Available with basket or shopping bag
- Maximum basket and shopping bag load 5kg / 11lbs
- 6 inch EVA wheel

Highlight Features



Large padded seat



Maximum load for shopping bag is 5kg/11lbs



Effortless braking system



RG70

RG70	Handle height (min - max)	Seat height	Length	Width	Width between handle	Net weight	Weight capacity
STANDARD	87.5 - 99.5cm / 34.4" - 39.2"	61cm / 24"	65cm / 25.6"	65cm / 25.6"	48cm / 18.9"	7.3kg / 16.1lbs	113kg / 250lbs
LOW	79.5 - 88.5cm / 31.3" - 34.8"	54cm / 21.3"	65cm / 25.6"	65cm / 25.6"	48cm / 18.9"	7.2kg / 15.8lbs	113kg / 250lbs
SUPER LOW	76.5 - 85.5cm / 30" - 33.7"	48cm / 18.9"	65cm / 25.6"	65cm / 25.6"	48cm / 18.9"	7.1kg / 15.6lbs	113kg / 250lbs

- 4 handle height adjustments.
- Available with basket or shopping bag.
- Maximum basket and shopping bag load 5kg / 11lbs.
- 200mm EVA wheel.

Highlight Features



Large padded seat



Simple handle adjustment



Inverse handle with effortless braking system.



RG72

RG72	Handle height (min - max)	Seat height	Length	Width	Width between handle	Net weight	Weight capacity
STANDARD	85.5 - 94.5cm / 33.7" - 37.2"	61cm / 24"	65cm / 25.6"	65cm / 25.6"	48cm / 18.9"	6.9kg / 15.2lbs	113kg / 250lbs
LOW	77.5 - 86.5cm / 30.5" - 34.1"	54cm / 21.3"	65cm / 25.6"	65cm / 25.6"	48cm / 18.9"	6.8kg / 15.0lbs	113kg / 250lbs
SUPER LOW	74 - 83cm / 29" - 32.7"	48cm / 18.9"	65cm / 25.6"	65cm / 25.6"	48cm / 18.9"	6.7kg / 14.7lbs	113kg / 250lbs

- 4 handle height adjustments.
- Available with basket or shopping bag.
- Maximum basket and shopping bag load 5kg / 11lbs.
- 200mm EVA wheel.

Highlight Features



Large padded seat



Uses star handle knob for handle height adjustment.



Inverse handle with ergonomic loop brake

NEW



RG77

RG77	Handle height (min - max)	Seat height	Length	Width	Width between handle	Net weight	Weight capacity
STANDARD	87.5 - 99.5cm / 34.4" - 39.2"	61cm / 24"	65cm / 25.6"	65cm / 25.6"	48cm / 18.9"	6.3kg / 13.6lbs	113kg / 250lbs
LOW	79.5 - 88.5cm / 31.3" - 34.8"	54cm / 21.3"	65cm / 25.6"	65cm / 25.6"	48cm / 18.9"	6.2kg / 13.4lbs	113kg / 250lbs
SUPER LOW	76.5 - 85.5cm / 30" - 33.7"	48cm / 18.9"	65cm / 25.6"	65cm / 25.6"	48cm / 18.9"	6.1kg / 13.2lbs	113kg / 250lbs

- 4 handle height adjustments.
- Available with shopping bag only.
- Maximum shopping bag load 5kg / 11lbs.
- 200mm EVA wheel.

Highlight Features



Net seat



Secondary frame lock



Folds compactly for storage. 25cm depth when folded