

ACCESSORIES / SPARE PART

We are also provide accessories for better life and spare part for your future replacement.



Accessories

Cane Holder :

As accessories we have a cane holder that is .placed on the side of handle. This allows the user to hang a cane on rollator and bring together.

Crutch Holder :

Also for accessories together with cane holder for support the cane

Slowdown brake :

Slow-down brakes allow you to control the rollator's speed on an incline. Afraid your walker might run away from you on a hillside?

Slow-down brakes let you set the brake tension, so the wheels roll in a controlled fashion.



Backstrap



Basket



Shopping bag



Cup / bottle holder



Cane holder



Crutch holder



Tray Holder



Slowdown brake